

**Rabbit Anti-Donkey IgG Heavy & Light Chain antibody {Biotin} (STJ99424)**

STJ99424

**GENERAL INFORMATION**

<b>Product Type</b>	Secondary antibodies
<b>Short Description</b>	Biotin-conjugated rabbit polyclonal anti-Donkey IgG Heavy & Light Chain secondary antibody. For use in most research applications.
<b>Applications</b>	WB/IHC/IF/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Donkey

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	2 mg/mL
<b>Conjugation</b>	Biotin
<b>Purification</b>	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
<b>Dilution Range</b>	WB 1:1000-10000
	ELISA 1:10000-20000
	IHC 1:50-100
	IF 1:50-1000
<b>Formulation</b>	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen Region</b>
<b>Specificity</b>
<b>Immunogen</b>
<b>Sequence</b>