

Goat Anti-Tree Shrew IgG Heavy & Light Chain antibody {Biotin} (STJ99438)

STJ99438

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Biotin-conjugated goat polyclonal anti-Tree Shrew IgG Heavy & Light Chain secondary antibody. For use in most research applications.
Applications	WB/IHC/IF/ELISA
Host/Source	Goat
Reactivity	Tree Shrew

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	2 mg/mL
Conjugation	Biotin
Purification	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Dilution Range	WB 1:1000-10000
	ELISA 1:10000-20000
	IHC 1:50-100
	IF 1:50-1000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence