

Goat Anti-Mouse IgG2a antibody {FITC} (STJ99449)

STJ99449

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	FITC-conjugated goat polyclonal anti-Mouse IgG2a secondary antibody. For use in most research applications.
Applications	WB/IHC/IF/ELISA
Host/Source	Goat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	2 mg/mL
Conjugation	FITC
Purification	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Dilution Range	IF 1:250-1:2000 FCM 1:250-1:2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence