

Goat Anti-Mouse IgG2b antibody {FITC} (STJ99450)

STJ99450

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Secondary antibodies |
| Short Description | FITC-conjugated goat polyclonal anti-Mouse IgG2b secondary antibody. For use in most research applications. |
| Applications | WB/IHC/IF/ELISA |
| Host/Source | Goat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 2 mg/mL |
| Conjugation | FITC |
| Purification | The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen. |
| Dilution Range | IF 1:250-1:2000 FCM 1:250-1:2000 |
| Formulation | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide. |
| Isotype | IgG2b |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence