

Rabbit Anti-Monkey IgG Heavy & Light Chain antibody {FITC} (STJ99463)

STJ99463

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	FITC-conjugated rabbit polyclonal anti-Monkey IgG Heavy & Light Chain secondary antibody. For use in most research applications.
Applications	WB/IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	2 mg/mL
Conjugation	FITC
Purification	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Dilution Range	IF 1:250-1:2000 FCM 1:250-1:2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence