

Rabbit Anti-Goat IgG Heavy & Light Chain antibody {Unconjugated} (STJ99567)

STJ99567

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Unconjugated rabbit polyclonal anti-Goat IgG Heavy & Light Chain secondary antibody. For use in most research applications.
Applications	WB/IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Goat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	2 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Dilution Range	WB 1:1000-10000 ELISA 1:10000-20000 IHC 1:50-100 IF 1:50-1000
Formulation	Liquid in PBS containing 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	