

Rabbit Anti-Guinea Pig IgG Heavy & Light Chain antibody {Unconjugated} (STJ99572)

STJ99572

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Unconjugated rabbit polyclonal anti-Guinea Pig IgG Heavy & Light Chain secondary antibody. For use in most research applications.
Applications	WB/IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Guinea Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	2 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Dilution Range	WB 1:1000-10000 ELISA 1:10000-20000 IHC 1:50-100 IF 1:50-1000
Formulation	Liquid in PBS containing 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence