

## Rabbit Anti-Duck IgY Heavy & Light Chain antibody {Unconjugated} (STJ99574)

STJ99574

### GENERAL INFORMATION

<b>Product Type</b>	Secondary antibodies
<b>Short Description</b>	Unconjugated rabbit polyclonal anti-Duck IgY Heavy & Light Chain secondary antibody. For use in most research applications.
<b>Applications</b>	WB/IHC/IF/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Duck

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	2 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
<b>Dilution Range</b>	WB 1:1000-10000 ELISA 1:10000-20000 IHC 1:50-100 IF 1:50-1000
<b>Formulation</b>	Liquid in PBS containing 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen Region</b>
<b>Specificity</b>
<b>Immunogen</b>
<b>Sequence</b>