

Anti-CD11c/Itgax antibody [N418] {PerCP/Cyanine5.5} (STJA0000068)

STJA0000068

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-Integrin alpha-X is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	N418
Concentration	5 Mu L/Test
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16411
Gene Symbol	Itgax
Uniprot ID	ITAX_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



C57BL/6 murine splenocytes were stained with Anti-Mouse CD11c Monoclonal Antibody (PerCP-Cy5.5 Conjugated) and Anti-Mouse MHC II (I-A/I-E) Monoclonal Antibody (FITC Conjugated) (left). Splenocytes stained with Anti-Mouse MHC II (I-A/I-E) Monoclonal Antibody (FITC Conjugated) (right) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081