

Anti-CD11c/Itgax antibody [N418] {APC} (STJA0000073)

STJA0000073

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-Integrin alpha-X is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	N418
Concentration	0.2 mg/mL
Conjugation	APC
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16411
Gene Symbol	Itgax
Uniprot ID	ITAX_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



C57BL/6 murine splenocytes were stained with Anti-CD11c antibody (STJA0000073) [Used at 0.2 μ g/ 10^6 cells dilution] and Anti-Mouse MHC II (I-A/I-E) Monoclonal Antibody (FITC Conjugated). Splenocytes stained with Anti-Mouse MHC II (I-A/I-E) Monoclonal Antibody (FITC Conjugated) (Right) were used as control

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081