

Anti-Cd80 antibody [16-10A1] {PerCP/Cyanine5.5} (STJA0000085)

STJA0000085

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-T-lymphocyte activation antigen CD80 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	16-10A1
Concentration	5 Mu L/Test
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12519
Gene Symbol	Cd80
Uniprot ID	CD80_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	



Mouse abdominal macrophages elicited by starch broth were stained with Anti-CD80 antibody (STJA0000085) (filled gray histogram). Unstained macrophages (blank black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081