

Anti-Adgre1 antibody [CI:A3-1] {AF488} (STJA0000129) STJA0000129

GENERAL INFORMATION

Product Type Primary antibodies Short Description Rat monoclonal antibody anti-Adhesion G protein-coupled receptor E1 is suitable for use in Flow Cytometry research applications. Applications FC Host/Source Rat

Reactivity Mouse

PRODUCT PROPERTIES


Clonality Monoclonal Clone ID CI:A3-1 Concentration 0.5 mg/mL Conjugation AF488 Purification
 Full induction

 Dilution Range

 Formulation
 Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Storage Instruction Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID 13733 Gene Symbol Adgre1 Uniprot ID AGRE1_MOUSE Immunogen Immunogen Region Specificity Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes. St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081