

Anti-CD3E antibody [OKT3] {PerCP} (STJA0000154)

STJA0000154

GENERAL INFORMATION


| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-T-cell surface glycoprotein CD3 epsilon chain is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | OKT3 |
| Concentration | 5 Mu L/Test |
| Conjugation | PerCP |
| Purification | |
| Dilution Range | |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG2ak |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|-------------------------|------------|
| Gene ID | 916 |
| Gene Symbol | CD3E |
| Uniprot ID | CD3E_HUMAN |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen | |
| Sequence | |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081