

Anti-CD3 antibody [17A2] {PE/Cyanine5} (STJA0000185)

STJA0000185

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD3 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse


PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	17A2
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID [12500](#)
[12502](#)
[12503](#)
[Cd3d](#)
[Cd3g](#)
[CD3D_MOUSE](#)
[CD3G_MOUSE](#)
[CD3Z_MOUSE](#)

Immunogen
Immunogen Region
Specificity
Immunogen Sequence



C57BL/6 murine splenocytes were stained with Anti-CD3 antibody (STJA0000185) (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.