

Anti-CD3 antibody [17A2] {PerCP/Cyanine5.5} (STJA0000187)

STJA0000187

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD3 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse


PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	17A2
Concentration	5 µg L/Test
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID [12502](#)
[12500](#)
[12503](#)
[Cd3g](#)
[Cd3d](#)
[CD3G_MOUSE](#)
[CD3D_MOUSE](#)
[CD3Z_MOUSE](#)

Immunogen
Immunogen Region
Specificity
Immunogen Sequence



C57BL/6 murine splenocytes were stained with Anti-CD3 antibody (STJA0000187) (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.