

Anti-CD3 antibody [17A2] {AF488} (STJA0000197)

STJA0000197

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD3 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse


PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	17A2
Concentration	0.5 mg/mL
Conjugation	AF488
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID [12502](#)
[12503](#)
[12501](#)
[Cd3g](#)
[Cd247](#)
[CD3G_MOUSE](#)
[CD3Z_MOUSE](#)
[CD3E_MOUSE](#)

Immunogen
Immunogen Region
Specificity
Immunogen Sequence



Mouse splenocytes were stained with Anti-CD3 antibody (STJA0000197) [Used at 0.2 Aug/10⁶ cells dilution] (filled gray histogram). Unstained splenocytes (blank black histogram) were used as control.