

Anti-CD3 antibody [17A2] {AF647} (STJA0000198)

STJA0000198

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD3 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse


PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	17A2
Concentration	0.5 mg/mL
Conjugation	AF647
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID [12501](#)
[12503](#)
[12500](#)
[Cd3e](#)
[Cd247](#)
[CD3E_MOUSE](#)
[CD3Z_MOUSE](#)
[CD3D_MOUSE](#)

Immunogen
Immunogen Region
Specificity
Immunogen Sequence



C57BL/6 murine splenocytes were stained with Anti-Mouse CD3 Monoclonal Antibody (AF647 Conjugated) (Used at 0.2 μ g/ 10^6 cells dilution) (filled gray histogram). Unstained lymphocytes (empty black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081