

Anti-CD47 antibody [MIAP410] {PerCP} (STJA0000210)

STJA0000210

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD47 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human/Mouse/Rat


PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MIAP410
Concentration	0.2 mg/mL
Conjugation	PerCP
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID [961](#)
[16423](#)
[29364](#)
[CD47](#)
[Cd47](#)
[CD47_HUMAN](#)
[CD47_MOUSE](#)
[CD47_RAT](#)

Immunogen
Immunogen Region
Specificity
Immunogen Sequence



C57BL/6 murine splenocytes were stained with Anti-CD47 antibody (STJA0000210) [Used at 0.2 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081