

Anti-Cd28 antibody [37.51] {PE} (STJA0000271)

STJA0000271

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Hamster monoclonal antibody anti-T-cell-specific surface glycoprotein CD28 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	37.51
Concentration	0.2 mg/mL
Conjugation	PE
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12487
Gene Symbol	Cd28
Uniprot ID	CD28_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



C57BL/6 murine splenocytes were stained with Anti-CD28 antibody (STJA0000271) [Used at 0.02 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081