

## Anti-Cd1d1 antibody [19G11] {FITC} (STJA0000295)

STJA0000295

### GENERAL INFORMATION


|                          |                                                                                                                                 |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                              |
| <b>Short Description</b> | Rat monoclonal antibody anti-antigen-presenting glycoprotein CD1d1 is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC                                                                                                                              |
| <b>Host/Source</b>       | Rat                                                                                                                             |
| <b>Reactivity</b>        | Mouse                                                                                                                           |

### PRODUCT PROPERTIES

|                            |                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Monoclonal                                                                                                   |
| <b>Clone ID</b>            | 19G11                                                                                                        |
| <b>Concentration</b>       | 5 Mu L/Test                                                                                                  |
| <b>Conjugation</b>         | FITC                                                                                                         |
| <b>Purification</b>        |                                                                                                              |
| <b>Dilution Range</b>      |                                                                                                              |
| <b>Formulation</b>         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.                  |
| <b>Isotype</b>             | IgG2bk                                                                                                       |
| <b>Storage Instruction</b> | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|                           |             |
|---------------------------|-------------|
| <b>Gene ID</b>            | 12479       |
| <b>Gene Symbol</b>        | Cd1d1       |
| <b>Uniprot ID</b>         | CD1D1_MOUSE |
| <b>Immunogen</b>          |             |
| <b>Immunogen Region</b>   |             |
| <b>Specificity</b>        |             |
| <b>Immunogen Sequence</b> |             |



C57BL/6 murine splenocytes were stained with Anti-CD1d1 antibody (STJA0000295) (filled gray histogram). Unstained lymphocytes (empty black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081