

Anti-CD11b antibody [M1/70] {AF488} (STJA0000437)

STJA0000437

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD11b is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M1/70
Concentration	0.5 mg/mL
Conjugation	AF488
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3684
	NA
Gene Symbol	ITGAM
	Itgam
Uniprot ID	ITAM_HUMAN
	ITAM_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Mouse bone marrow cells were stained with Anti-Mouse/Human CD11b antibody (STJA0000437) [Used at 0.2 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained bone marrow cells (blank black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081