

Anti-CD25/Il2ra antibody [PC-61.5.3] {APC} (STJA0000570)

STJA0000570

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Interleukin-2 receptor subunit alpha is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PC-61.5.3
Concentration	0.2 mg/mL
Conjugation	APC
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16184
Gene Symbol	Il2ra
Uniprot ID	IL2RA_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



C57BL/6 murine splenocytes were stained with Anti-CD25 antibody (STJA0000570) [Used at 0.2 μ g/ 10^6 cells dilution] and Anti-Mouse CD4 Monoclonal Antibody (FITC Conjugated). Splenocytes stained with Anti-Mouse CD4 Monoclonal Antibody (FITC Conjugated) and Rat IgG1 Isotype Control (APC Conjugated) (right) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081