

Anti-CD45R/B220 antibody [RA3.3A 1/6.1] {PerCP/Cyanine5.5} (STJA0000674)

STJA0000674

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD45R/B220 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RA3.3A 1/6.1
Concentration	0.2 mg/mL
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgMk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	5788
	19264
Gene Symbol	PTPRC
	Ptprc
Uniprot ID	PTPRC_HUMAN
	PTPRC_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



C57BL/6 murine splenocytes were stained with Anti-CD45R/B220 antibody (STJA0000674) [Used at 0.2 µg/10⁶ cells dilution] (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081