

Anti-CD49b/Itga2 antibody [DX5] {AF488} (STJA0000703)

STJA0000703

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Integrin alpha-2 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DX5
Concentration	0.5 mg/mL
Conjugation	AF488
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgMk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16398
Gene Symbol	Itga2
Uniprot ID	ITA2_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Mouse bone marrow cells were stained with Anti-CD49b antibody (STJA0000703) [Used at 0.2 Aug/10⁶ cells dilution] (filled gray histogram). Unstained bone marrow cells (blank black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081