

Anti-TCR Beta antibody [H57-597 (HB218)] {AF647} (STJA0000783)

STJA0000783

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-TCR Beta is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	H57-597 (HB218)
Concentration	0.5 mg/mL
Conjugation	AF647
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081