

Anti-TER-119 antibody [TER-119] {PerCP/Cyanine5.5} (STJA0000788)

STJA0000788

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-TER-119 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TER-119
Concentration	5 Mu L/Test
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence



C57BL/6 murine bone marrow cells were stained with Anti-TER-119 antibody (STJA0000788) (filled gray histogram). Unstained bone marrow cells (empty black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081