

Anti-CD34 antibody [581] {PE} (STJA0000903)

STJA0000903

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Hematopoietic progenitor cell antigen CD34 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	581
Concentration	5 Mu L/Test
Conjugation	PE
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	947
Gene Symbol	CD34
Uniprot ID	CD34_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Human T lymphoblastic leukemia cell line Molt-4 were stained with Anti-CD34 antibody (STJA0000903) (filled gray histogram) or Mouse IgG1 isotype Control PE (empty black histogram).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081