

## Anti-Cd14 antibody [Sa14-2] {PE} (STJA0000979)

STJA0000979

### GENERAL INFORMATION


|                          |                                                                                                                                 |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                              |
| <b>Short Description</b> | Rat monoclonal antibody anti-Monocyte differentiation antigen CD14 is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC                                                                                                                              |
| <b>Host/Source</b>       | Rat                                                                                                                             |
| <b>Reactivity</b>        | Mouse                                                                                                                           |

### PRODUCT PROPERTIES

|                            |                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Monoclonal                                                                                                   |
| <b>Clone ID</b>            | Sa14-2                                                                                                       |
| <b>Concentration</b>       | 5 Mu L/Test                                                                                                  |
| <b>Conjugation</b>         | PE                                                                                                           |
| <b>Purification</b>        |                                                                                                              |
| <b>Dilution Range</b>      |                                                                                                              |
| <b>Formulation</b>         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.                  |
| <b>Isotype</b>             | IgG2ak                                                                                                       |
| <b>Storage Instruction</b> | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|                           |            |
|---------------------------|------------|
| <b>Gene ID</b>            | 12475      |
| <b>Gene Symbol</b>        | Cd14       |
| <b>Uniprot ID</b>         | CD14_MOUSE |
| <b>Immunogen</b>          |            |
| <b>Immunogen Region</b>   |            |
| <b>Specificity</b>        |            |
| <b>Immunogen Sequence</b> |            |



Mouse abdominal macrophages elicited by starch broth were stained with Anti-CD14 antibody (STJA0000979) (filled gray histogram). Unstained macrophages (blank black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081