

Anti-CD90.1 antibody [OX-7] {PerCP/Cyanine5.5} (STJA0001084)

STJA0001084

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD90/CD90.1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OX-7
Concentration	0.2 mg/mL
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	24832
Gene Symbol	Thy1
Uniprot ID	THY1_RAT
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Rat thymocytes were stained with Anti-CD90/Mouse CD90.1 antibody (STJA0001084) [Used at 0.2 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained thymocytes (empty black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081