

Anti-Cd200 antibody [OX-90] {APC} (STJA0001116)

STJA0001116

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-OX-2 membrane glycoprotein is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OX-90
Concentration	0.2 mg/mL
Conjugation	APC
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	
Gene Symbol	Cd200
Uniprot ID	OX2G_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



C57BL/6 murine splenocytes were stained with Anti-CD200 antibody (STJA0001116) and Anti-Mouse CD45R/B220 Monoclonal Antibody (FITC Conjugated).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081