

Anti-Cd19 antibody [1D3] (STJA0001210)

STJA0001210

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-B-lymphocyte antigen CD19 is suitable for use in Universal, Universal, Universal, Flow Cytometry and Immunoprecipitation research applications.
Applications	Cell Sep-Neg/Depletion/Neut/FC/IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1D3
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2ak
Storage	Store at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	12478
Gene Symbol	Cd19
Uniprot ID	CD19_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	