

Anti-Cd80 antibody [16-10A1] (STJA0001228)

STJA0001228

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-T-lymphocyte activation antigen CD80 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
Applications	FC/IHC/IP
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	16-10A1
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG
Storage	Store at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	12519
Gene Symbol	Cd80
Uniprot ID	CD80_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	