

Anti-CD3E antibody [OKT3] (STJA0001257)

STJA0001257

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-T-cell surface glycoprotein CD3 epsilon chain is suitable for use in Flow Cytometry and Immunohistochemistry research applications.
Applications	FC/IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OKT3
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2ak
Storage	Store at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	916
Gene Symbol	CD3E
Uniprot ID	CD3E_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	