

Anti-CD3 antibody [17A2] (STJA0001276)

STJA0001276

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD3 is suitable for use in Flow Cytometry, Immunohistochemistry, Immunoprecipitation and Immunocytochemistry research applications.
Applications	FC/IHC-F/IP/ICC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	17A2
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2bk
Storage	Store at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	12503 12500 12501 Cd247 Cd3d CD3Z_MOUSE CD3D_MOUSE CD3E_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	