

Anti-CD3 antibody [17A2] (STJA0001277)

STJA0001277

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD3 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
Applications	FC/IHC/IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	17A2
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2bk
Storage Instruction	Store at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	12501 12503 12502 Cd3e Cd247 CD3E_MOUSE CD3Z_MOUSE CD3G_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	