

Anti-CD132/Il2rg antibody [3E12] {Biotin} (STJA0001283)

STJA0001283

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Cytokine receptor common subunit gamma is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3E12
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	16186
Gene Symbol	IL2rg
Uniprot ID	IL2RG_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	