

Anti-CD47 antibody [MIAP410] (STJA0001288)

STJA0001288

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD47 is suitable for use in Flow Cytometry, Western Blot, Immunohistochemistry, Immunoprecipitation, Immunocytochemistry and Immunofluorescence research applications.
Applications	FC/WB/IHC/IP/ICC/IF
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MIAP410
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG1k
Storage	Store at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	16423 961 29364 Cd47 CD47 CD47_MOUSE CD47_HUMAN CD47_RAT
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	