

Anti-CD47 antibody [MIAP410] {PE/Cyanine5.5} (STJA0001292)

STJA0001292

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD47 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MIAP410
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	29364 961 16423 Cd47 CD47 CD47_RAT CD47_HUMAN CD47_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	