

Anti-Cd48 antibody [HM48-1] (STJA0001295)

STJA0001295

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-CD48 antigen is suitable for use in Universal, Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
Applications	Block/FC/IHC/IP
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HM48-1
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG
Storage	Store at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	12506
Gene Symbol	Cd48
Uniprot ID	CD48_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	