

**Anti-Cd48 antibody [HM48-1] {Biotin} (STJA0001297)**

STJA0001297

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Armenian Hamster monoclonal antibody anti-CD48 antigen is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Armenian Hamster
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	HM48-1
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Biotin
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at 2-8°C. Do not freeze this product.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">12506</a>
<b>Gene Symbol</b>	<a href="#">Cd48</a>
<b>Uniprot ID</b>	<a href="#">CD48_MOUSE</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	