

**Anti-Cd22 antibody [Cy34.1] (STJA0001314)**

STJA0001314

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-B-cell receptor CD22 is suitable for use in Flow Cytometry and Immunoprecipitation research applications.
<b>Applications</b>	FC/IP
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	Cy34.1
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.05% Proclin300
<b>Isotype</b>	IgG1k
<b>Storage Instruction</b>	Store at 2-8°C. Do not freeze this product.

**TARGET INFORMATION**

<b>Gene ID</b>	12483
<b>Gene Symbol</b>	Cd22
<b>Uniprot ID</b>	CD22_MOUSE
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	