

Anti-Cd22 antibody [Cy34.1] {PE/Cyanine5.5} (STJA0001326)
STJA0001326**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-B-cell receptor CD22 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Cy34.1
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12483
Gene Symbol	Cd22
Uniprot ID	CD22_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	