

**Anti-Cd40 antibody [FGK4.5/FGK45] {Biotin} (STJA0001354)**

STJA0001354

**GENERAL INFORMATION**

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Rat monoclonal antibody anti-Tumor necrosis factor receptor superfamily member 5 is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC  |
| <b>Host/Source</b>       | Rat   |
| <b>Reactivity</b>        | Mouse   |

**PRODUCT PROPERTIES**

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Monoclonal  |
| <b>Clone ID</b>            | FGK4.5/FGK45  |
| <b>Concentration</b>       | 0.5 mg/mL   |
| <b>Conjugation</b>         | Biotin  |
| <b>Purification</b>        |   |
| <b>Dilution Range</b>      |   |
| <b>Formulation</b>         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| <b>Isotype</b>             | IgG2ak  |
| <b>Storage Instruction</b> | Store at 2-8°C. Do not freeze this product.   |

**TARGET INFORMATION**

|                    |                             |
|--------------------|-----------------------------|
| <b>Gene ID</b>     | <a href="#">21939</a>       |
| <b>Gene Symbol</b> | <a href="#">Cd40</a>        |
| <b>Uniprot ID</b>  | <a href="#">TNFR5_MOUSE</a> |
| <b>Immunogen</b>   |                             |
| <b>Immunogen</b>   |                             |
| <b>Region</b>      |                             |
| <b>Specificity</b> |                             |
| <b>Immunogen</b>   |                             |
| <b>Sequence</b>    |                             |