

Anti-Cd40 antibody [FGK4.5/FGK45] {PE/Cyanine5.5} (STJA0001357)

STJA0001357

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Tumor necrosis factor receptor superfamily member 5 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	FGK4.5/FGK45
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	21939
Gene Symbol	Cd40
Uniprot ID	TNFR5_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	