

Anti-CD122/IL2rb antibody [5H4] (PE) (STJA0001372)

STJA0001372

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Interleukin-2 receptor subunit beta is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5H4
Concentration	0.2 mg/mL
Conjugation	PE
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16185
Gene Symbol	IL2rb
Uniprot ID	IL2RB_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	