

Anti-Cd1d1 antibody [19G11] {PE/Cyanine5.5} (STJA0001386)

STJA0001386

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-antigen-presenting glycoprotein CD1d1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	19G11
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12479
Gene Symbol	Cd1d1
Uniprot ID	CD1D1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	