

Anti-CD29/ITGB1 antibody [TS2/16.2.1] (STJA0001434)

STJA0001434

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Integrin beta-1 is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Immunoprecipitation research applications. |
| Applications | FC/IHC/IF/IP |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-----------------------|---|
| Clonality | Monoclonal |
| Clone ID | TS2/16.2.1 |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300 |
| Isotype | IgG1k |
| Storage | Store at 2-8°C. Do not freeze this product. |
| Instruction | |

TARGET INFORMATION

| | |
|--------------------|----------------------------|
| Gene ID | 3688 |
| Gene Symbol | ITGB1 |
| Uniprot ID | ITB1_HUMAN |
| Immunogen | |
| Immunogen | |
| Region | |
| Specificity | |
| Immunogen | |
| Sequence | |