

Anti-CD18/ITGB2 antibody [TS1/18.1.2.11] {PE/Cyanine5} (STJA0001464)

STJA0001464

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Integrin beta-2 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TS1/18.1.2.11
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3689
Gene Symbol	ITGB2
Uniprot ID	ITB2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	