

Anti-CD57/B3GAT1 antibody [HNK-1] {Biotin} (STJA0001486)

STJA0001486

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Galactosylgalactosylxylosylprotein 3-beta-glucuronosyltransferase 1 is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-----------------------|---|
| Clonality | Monoclonal |
| Clone ID | HNK-1 |
| Concentration | 0.5 mg/mL |
| Conjugation | Biotin |
| Purification | |
| Dilution Range | |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgMk |
| Storage | Store at 2-8°C. Do not freeze this product. |
| Instruction | |

TARGET INFORMATION

| | |
|--------------------|-----------------------------|
| Gene ID | 27087 |
| Gene Symbol | B3GAT1 |
| Uniprot ID | B3GA1_HUMAN |
| Immunogen | |
| Immunogen | |
| Region | |
| Specificity | |
| Immunogen | |
| Sequence | |