

Anti-CD32/FCGR2B antibody [IV-3] {PE/Cyanine7} (STJA0001505)

STJA0001505

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Low affinity immunoglobulin gamma Fc region receptor II-b is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	IV-3
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine7
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	2213
Gene Symbol	FCGR2B
Uniprot ID	FCG2B_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	